

## HDL - Cholesterol

### Precept. Reagent


HDL-MC-0150 (1x50ml) HDL-MC-0520 (5x20ml)  
HDL-MC-0620 (6x20ml)

### INTENDED FOR USE:

For the quantitative determination of HDL cholesterol in serum and plasma

### PRINCIPLE :

HDL-Cholesterol is obtained through selective precipitation of LDL and VLDL lipoproteins, thus HDL lipoproteins remain in solution. HDL-Cholesterol in supernatant is treated as a sample for cholesterol assay according to the following reaction:



Formed color is measured at 546 nm and is proportional to HDL-Cholesterol concentration in sample when used as directed.

### SPECIMEN COLLECTION:

Fresh ( or just defrosted ) not hemolyzed serum or plasma ( EDTA Na<sub>2</sub> , Na Heparin ).

Centrifuge and collect serum as soon as possible .

HDL<sub>C</sub> in serum or plasma is stable up to 7 days at +2-8°C ,1 month at 20°C and 2 years at -70°C Shake and bring the samples at room temperature ( +15-25°C. ) before using .

### REAGENTS COMPOSITION :

The reagents set are stored at ambient temperature. Storage must not exceed expiration date on box label.

Precipitating Reagent	
Phosphotungstic acid	0.55 mM
Magnesium Chloride	25 mM

### PACKAGE : Collection & Storage .

Store at +2-8°C.

Stable until the expiration date reported upon the package.

After the unsealing and the taking of the reagent , it is advised to close up the bottle immediately in order to avoid evaporation , direct light exposure and bacterial contamination .

### PRECAUTION & WARNING:

Avoid pipetting by mouth .

The preparation , according to current regulation . is classified as not dangerous.

The total concentration of non active components ( preservatives , detergents ,stabilizers ) is below the minimum required for citation .

Anyway handle with care , avoid ingestion , avoid contact with eyes , skin and mucous membranes . The samples must be handle as potentially infected from HIV or Hepatitis .

### REQUIRED MATERIALS NOT PROVIDED :

General Laboratory Equipment and instrumentations .

### PROCEDURE :

This methodology describes the manual procedure to use the kit .

For automated procedure, ask for specific application .

### Precipitation

Specimen	200 µl
Precipitant	500 µl

Mix and allow standing for 10 minutes at room temperature. Centrifuge for 10 minutes at 4000 rpm, or 2 minutes 12000rpm. Separate of the clear supernatant within two hours and determine the cholesterol content by the CHOD-PAP method. The supernatant may be stored up to five days at 2-8°C

Wavelength 546 nm

Optical path 1 cm

Incubation temperature 20, 25 or 37°C

Zero adjustment Reagent blank

	Blank	Standard	Specimen
cholesterol reagent	1.0 ml	1.0 ml	1.0 ml
Specimen	—	—	100 µl
Distilled water	100 µl	—	—
Standard (R1)	—	100 µl	—
	Mix, and incubate for 5 min. at 37°C or 10 min. at 20-25°C. Measure absorbance of specimen (A <sub>specimen</sub> ) against reagent blank within 30 minutes.		

### CALCULATION :

#### 1. HDL Cholesterol

$$\text{mg/dl} = 50 \times \frac{(\text{A}) \text{ Sample}}{(\text{A}) \text{ Standard}}$$

#### 2. LDL Cholesterol

$$\text{LDL Cholesterol (mg/dl)} = \text{Total Cholesterol} - \frac{\text{Triglycerides}}{5} \text{ HDL Cholesterol}$$

### EXPECTED VALUES :

<b>HDL:</b>	
<b>Women</b>	30-85 mg/dL
<b>Men</b>	30-70 mg/dL
<b>LDL: Adults</b>	66-178 mg/dL

The above mentioned values are to be considered as a reference.

It is strongly recommended that each laboratory establish its own normal range according to its geographic area , according to IFCC protocol .

Since CHOL HDL has an elevated protective action against the risk of arising cardiovascular diseases , the following reference values can be used :

Protective Action	Men	Women
High	> 55mg/dL	> 65 mg/dL
None	35-55 mg/dL	45-65 mg/dL
Poor	< 35 mg/dL	< 45 mg/dL

### WASTE DISPOSAL :

The disposal of the product must be in accordance with local regulation concerning waste disposal

### QUALITY CONTROL :

It is recommended to execute the quality control at every kit utilization to verify that values are within the reference range indicated by the methodology.

### PERFORMANCE :

MEASURE INTERVAL LINEARITY :	2-200 mg/dl
LOWEST MEASURABLE LIMIT :	2 mg/dl
SENSITIVITY :	1mg/dl= ΔA

### PRECISION WITH SERIES : n=20

LEVEL	M	C.V
LOW LEVEL	M = 30.2 mg/dl	C.V = 1.50%
MEDIUM LEVEL	M = 42.7 mg/dl	C.V = 1.18%
HIGH LEVEL	M = 75.6 mg/dl	C.V = 1.04%

### PRECISION AMONG SERIES : n=20

LEVEL	M	C.V
LOW LEVEL	M = 29.4 mg/dl	C.V = 2.68%
MEDIUM LEVEL	M = 41.9 mg/dl	C.V = 1.89%
HIGH LEVEL	M = 73.2 mg/dl	C.V = 3.22%
INTER. ANALYZED	30-72 mg/dl	
CORRELATION	r = 0.999	n=50
LIN. REGRESSION	y= 1.01 x- 3.39	n=50

### INTERFERENCE:

Interferences are negligible up to :

Bilirubin	20 mg/dl		
Hemoglobin	0.4 g/dl	Glucose	500 mg/dl
Ascorbic Acid	40 mg/dl	Triglycerids	2000 mg/dl

### METHOD LIMITATIONS:

If Triglycerides levels are higher than 2000 mg/dl , repeat the measure on a sample diluted 1:2 with physiological solution e multiply the results × 2 .

Do not use Anticoagulants containing citrate .

For thorough evaluation of the interfering substances ,consult : Young , D. S .et al , AACC Press , Washington DC , , 3-104 ( 1990 )

### REFERENCES:

1. Tietz, N.W. (ed) Fundamentals of Clinical Chemistry W.B.Saunders Co., Philadelphia, 1976.
2. Watson, D., Clin. Chem. Acta 5 (637), 1960.
3. Trinder, P., Ann Clin. Biochem. 6 (24), 1969.
4. Castelli, W.P., et al., Circ. 55 (767) 1977.

 <p><b>NS BIOTEC</b> Medical EQUIPMENT</p> <p>66 port said st. camp shezar -Alexandria- Egypt Tel : 002 03 592 0902 Fax: 002 03 592 0908 Website: www.nsbiotec.com Email: info@nsbiotec.com</p>	 <p><b>CE</b> <b>EC REP</b></p> <p>.CMC Medical Devices &amp; Drugs S.L C/ Horacio Lengo, 18. 29006. Málaga, Spain</p>
--	---