

ALT/ALAT/GPT (4+1)

Kinetic

REF. GPT –MK – 0520 (5 X 20 ml)

. GPT –MK – 1010 (10X10 ml)

INTENDED USE

NS Biotec ALT reagent is intended for the in vitro quantitative determination of alanine aminotransferase (EC 2.6.1.2) activity in serum on both automated and manual systems.

CLINICAL SIGNIFICANCE

Alanine aminotransferase (glutamate pyruvate transaminase) belongs to the group of transaminases, which catalyze the conversion of amino acids to the corresponding α -keto acids via the transfer of amino groups; they also catalyze the reverse process. Although higher activities exist in the liver, minor activity can also be detected in the kidneys, heart, skeletal muscle, pancreas, spleen, and lungs. Elevated serum ALT is found in hepatitis, cirrhosis, obstructive jaundice, carcinoma of the liver, and chronic alcohol abuse. ALT is only slightly elevated in patients who have an uncomplicated myocardial infarction. Although both serum aspartate aminotransferase (AST) and ALT become elevated whenever disease processes affect liver cell integrity, ALT is the more liver-specific enzyme. Moreover, elevations of ALT persist longer than elevations of AST activity.

ASSAY PRINCIPLE

ALT methods using U.V. procedures were first described by Henley in 1955. Wroblewski and LaDue in 1956 described a method for determining ALT using lactate dehydrogenase (LDH) and nicotinamide adenine dinucleotide, reduced (NADH). This method was later modified by Henry and Bergmeyer to optimize substrate conditions and eliminate side reactions. The NS Biotec ALT reagent is based on the recommendations of the IFCC. The series of reactions involved in the assay system is as follows:

1. The amino group is enzymatically transferred by ALT present in the specimen from alanine to the carbon atom of 2-oxoglutarate yielding pyruvate and L-glutamate.
2. Pyruvate is reduced to lactate by LDH present in the reagent with the simultaneous oxidation of NADH to NAD.



The rate of oxidation of the coenzyme NADH is proportional to the ALT activity in the specimen. It is determined by measuring the decrease in absorbance at 334 / 340 / 365 nm correspondingly. Endogenous specimen pyruvate is rapidly and completely reduced by LDH during the initial incubation period so that it does not interfere with the assay.

EXPECTED VALUES

	Male	Female
25°C	Up to 22 U/l Up to 0.37 μ kat/l	Up to 17 U/l Up to 0.28 μ kat/l
30°C	Up to 29 U/l Up to 0.48 μ kat/l	Up to 22 U/l Up to 0.37 μ kat/l
37°C	Up to 41 U/l Up to 0.68 μ kat/l	Up to 31 U/l Up to 0.52 μ kat/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the ALT results should

always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

REAGENTS

	Tris buffer pH 7.8	100 mmol/l
R ₁	L-Alanine	500 mmol/l
	LDH	1700 U/l
R ₂	NADH	0.18 mmol/l
	2-Oxoglutarate	18 mmol/l

• Reagent Preparation & Stability

All reagents are ready for use and stable up to the expiry date given on label when stored at 2–8°C.

Prepare the working solution by as follows:

The working reagent can be prepared by mixing 4 volume of reagent 1 (R1) and 1 volume of reagent 2 (R2).

- Stability after preparation.

2 days	at 15–25 °C.
4 weeks	at 2 – 8 °C

SPECIMEN

Serum, EDTA or heparinized plasma. Avoid hemolysis.

• Specimen Preparation & Stability

Separate serum/plasma from clot/cells within 8 hours at room temperature or 48 hours at 2–8°C.

Freezing of the samples is not recommended.

PROCEDURE

• Manual Procedure

Wavelength	340,334,365 nm
Cuvette	1 cm light path
Temperature	25, 30 or 37 °C
Zero adjustment	against air

Pipette into test tube or cuvette	
Working solution	1000 μ l
Serum or plasma	100 μ l

Mix, incubate for 1.0 minute, and start stopwatch simultaneously. Read again after exactly 1, 2, and 3 minutes.

• Automated Procedure

User defined parameters for different auto analyzers are available upon request.

CALCULATION

Determine the change in absorbance per minute ($\Delta A/\text{min}$) from the linear portion of the reaction curve and calculate the ALT/GPT activity by using the following formulae:

$$\begin{aligned} \text{U/l} &= 1780 \times \Delta A \text{ 334 nm/min} \\ \text{U/l} &= 1746 \times \Delta A \text{ 340 nm/min} \\ \text{U/l} &= 3235 \times \Delta A \text{ 365 nm/min} \end{aligned}$$

One international unit (**U**) is defined as the amount of enzyme that catalyzes the transformation of one micromole of substrate per minute under specified conditions.

The general formula for converting $\Delta A/\text{min}$ into U/l is:

$$U/l = \frac{\Delta A/\text{min} \times TV \times 1000}{\sum * \times SV \times LP}$$

Where:

TV	Total reaction volume in ml
SV	Sample volume in ml
$\sum *$	millimolar absorptivity of NADH
LP	Cuvette pathlength in cm
1000	Conversion of U/ml to U/l

* millimolar absorptivity of NADH

at 334 nm = 6.18,

at 340 nm = 6.22, and

at 365 nm = 6.40

• Unit conversion

$U/l \times 16.67 \times 10^{-3} = \mu\text{kat/l}$

• Temperature correction

Multiply the result by 1.31 if the assay performed at 25°C but is to be reported at 30°C.

Multiply the result by 1.91 if the assay performed at 25°C but is to be reported at 37°C.

Multiply the result by 1.43 if the assay performed at 30°C but is to be reported at 37°C.

LINEARITY

When run as recommended, the assay is linear up to 400 U/l or 6.63 $\mu\text{kat/l}$

If result exceeds 400 U/l or 6.63 $\mu\text{kat/l}$, specimen should be diluted 1+5 with 0.9% NaCl solution and reassayed. Multiply the result by 6.

SENSITIVITY

The sensitivity is defined as the lower detection limit represents the lowest measurable ALT/GPT activity that can be distinguished from zero.

When run as recommended the sensitivity of this assay is 5 U/l or 0.75 $\mu\text{kat/l}$

QUALITY CONTROL

It is recommended that controls (normal and abnormal) be included in:

- Each set of assays, or
- At least once a shift, or
- When a new bottle of reagent is used, or
- After preventive maintenance is performed or a clinical component is replaced.

Commercially available control material with established ALT/GPT values may be routinely used for quality control.

Failure to obtain the proper range of values in the assay of control material may indicate:

- Reagent deterioration,
- Instrument malfunction, or
- Procedure errors.

The following corrective actions are recommended in such situations:

- Repeat the same controls.

- If repeated control results are outside the limits, prepare fresh control serum and repeat the test.
- If results on fresh control material still remain outside the limits, then repeat the test with fresh reagent.
- If results are still out of control, contact NS Biotec Technical Services.

INTERFERING SUBSTANCES

• Anticoagulants:

Fluoride and citrate inhibit the enzyme activity. The only accepted anticoagulants are heparin and EDTA.

• Bilirubin:

No interference from free bilirubin up to a level of 15 mg/dl, and from conjugated bilirubin up to level of 6.8 mg/dl.

• Drugs:

Young⁷ in 1990 has published a comprehensive list of drugs and substances which may interfere with this assay.

• Haemolysis:

Erythrocyte contamination may elevate results, since ALT/GPT activities in erythrocytes are three to five times higher than those in normal sera.

• Lipemia:

Lipemic specimens may cause high absorbance flagging. Choose diluted sample treatment for automatic rerun.

• Pyruvate:










High levels of serum pyruvate may interfere with assay performance.

WARNING & PRECAUTIONS

- NS Biotec ALT reagent is for in vitro diagnostic use only. Normal precautions exercised in handling laboratory reagents should be followed.
- Warm up working solution to the corresponding temperature before use.
- The reagent and sample volumes may be altered proportionally to accommodate different spectrophotometer requirements.
- Valid results depend on an accurately calibrated instrument, timing, and temperature control.
- Don't use the reagent if it is turbid or if the absorbance is less than 1.0 at 340 nm.

BIBLIOGRAPHY

1. **Zilva, JF, Pannall, PR (1979):** Plasma enzymes in Diagnosis in Clinical Chemistry in Diagnosis and treatment. Lloyd-Luke London. 1979: Chap 17:338.
2. **Henley, KS, Pollard, HM (1955):** J. Lab Clin. Med., Vol. 46, p. 785.
3. **Wroblewski, F, LaDue, JS (1956):** Proc Sec Exp Biol and Med 1956; 34:381.
4. **Henry, RJ, et al. (1960):** Am J Clin Path 1960; 34: 381.
5. **Bergmeyer, HU, et al. (1978):** Clin chem. 1978; 24: 58-73.
6. **IFCC Export panel on enzymes (1986):** part 3. J Clin Chem Clin Biochem 1986; 24:481-95.
7. **Young, Ds (1990):** Effects of Drugs on Clinical Laboratory Tests. Third Edition. 1990: 3: 6-12.

	Consult Instruction for Use
	Caution Consult Accompanying Documents
	In Vitro Diagnostic Medical Device
	Temperature Limitation
	Manufacturer
	Authorized Representative In The European Community
	Catalogue Number
	Batch Code
	Use By

 <p>N.S BIOTEC MEDICAL EQUIPMENTS 66 Port Said St., Camp Shezar Alexandria – Egypt Tele: 002 03 592 0902 Fax : 002 03 592 0908 Website : www.nsbiotec.com E- mail : info@nsbiotec.com</p>	  <p>CMC Medical Devices & Drugs S.L. C/ Horacio Lengo, 18. 29006. Málaga, Spain</p>
---	---