




# Vertical Laminar Air Flow

Model: LF-600



Custom - Made  
Design  
& OEM Available

**NS. BIOTEC**  
[www.nsbiotec.com](http://www.nsbiotec.com)



# Vertical Laminar Air Flow

## Model: LF - 600



- **External Made of Structured Stainless Steel**
- **Internal made of Stainless Steel 304 B2**
- **Air Flow : 0.2-0.6 m/s**
- **Air flow control: Yes**
- **HEPA Filter: 99.999 % Efficiency at 0.3 pm**
- **Noise: < 60 DB (A)**
- **Power supply: AC220v / 50-60HZ**
- **Power 600 W — 3 A**
- **Dimension of working area (W•D•H): 545 \*400'440 mm**
- **Overall dimension (W\*D\*H): 550 \*460'1000 mm**
- **Ultraviolet light: Wave Length 250-280 , Life time 11000 Hr**
- **Normal Light: Led White**
- **Net Weight: 35 kg**
- **Gross Weight : 40 KG**
- **Packing Dimensions: 560 • 470 • 1010 mm**



**ISO  
13485**

Some of the illustrations in this publication include the special accessories; we reserve the right to make changes in technical specification. Dimension without prior notice.

Headquarters : 66 Port Said St., Camp Shezar. Alexandria, Egypt  
Factory: 1 Elmadina st, Indstrial area, Elbetash, Alexandria, Egypt  
E-mail: [info@nsbiotec.com](mailto:info@nsbiotec.com) [export@nsbiotec.com](mailto:export@nsbiotec.com)

Tel : 002 03 5920902 Fax: 002 03 5920908

[www.nsbiootec.com](http://www.nsbiootec.com)