



Orbital Shaker OS - 250



Custom - Made
Design
& OEM Available

NS. BIOTECH
www.nsbiotec.com



Orbital Shaker

Model: OS - 250



Description

Nsbiotec Orbital shaker is new release with all outstanding features, it's made of Stainless Steel 304. Digital control, high Loading capability, different platform selection for many applications very competitive pricing, very low noise level. Speed is up to 250 rpm Orbital movement stroke length about 20 mm .Timer for hours, minutes and seconds

Features

- Continuous one-way orbital shaking One platform
- 50kg carry capability
- Microprocessor control
- Timer

Specification

- Motion: Orbital
- Max. Stroke Length: 20mm
- Speed / Inc: 0~250rpm / 1 rpm
- Timer: 1~9999 mins with alarm, continuous, 1 min increment
- Operating Temperature: Ambient to 40°C
- Carry Capability: 50kg
- Platform Dimension: 400 x 300 mm (W x L)
- Additional Platform: Yes optional
- Platform Material: 304 Stainless steel
- Unit Dimension: 450 x 450 x 290 mm (W x L x H)
- Weight: approx. 7.0kg
- Rated Voltages: 220V - 50 HZ



ISO
13485

Universal platform with adjustable bars with non - slip rubber mat Suitable for use in different laboratories & blood bank

