




Manual Staining MSS 12 / MSS 3



SET WITH 12 STAINING RESERVOIRS



Rack 24 Slides



**SET WITH 3
STAINING RESERVOIRS**

NS. BIOTEC
www.nsbiotec.com




SET WITH 12 STAINING RESERVOIRS

Structure entirely made from stainless steel sheet, resistant To corrosion of main solvents used in anatomic Pathology as xylol and xylene alternative , complete with slide holder for 24 slides

Code: MSS – 12

- Material : stainless steel AISI 304
- Steel surface finishing : scotch bright or structured stainless steel 304
- Dimensions (L X W X H) : 685 x 130 x 125 mm
- Weight (with reservoirs): 6.5 kg.
- Lid holder
- 3 available timer models (optional)
- Timer support (optional)




SET WITH 3 STAINING RESERVOIRS

Model with 3 staining reservoirs, with slide rack for 24 slides , stainless Steel structure AISI 304 .

Code: MSS – 3

- Material : stainless steel AISI 304
- Steel surface finishing : scotch bright or structured stainless steel
- Dimensions (L X W X H) : 200 x 130 x 125 mm
- Weight (with reservoirs): 2.0 kg.
- Lid holder
- 3 available timer models (optional)
- Timer support (optional)



STAINING RESERVOIRS FOR REAGENTS.

Both the manual staining sets mentioned above, use reservoirs Made from ABS, resistant to The most common solvents & alcohols Used in anatomic pathology Labs as xylol & xylene Alternative

Code : MSB-1

- Material: ABS
- Dimensions (L X W X H) 107 x 50 x 98 mm



RACK FOR SLIDES

ABS rack with handle, able to hold up to 24 slides, available for both manual staining sets mentioned above.

Code : MSR-24

- Material: ABS
- Dimensions (L X W X H) 88 x 31 x 90 mm
- Capacity: 24 slides

