

*Crucial where
zero TDS is a must*



ULTRA-PURE WATER SYSTEM

RO-100

97% - 99%
HIGH DESALINIZATION RATE



① 66 Port Said st., - Camp Chezar - Alexandria - Egypt
☎ +2 03 5920 902
🌐 www.nsbiotec.com

☎ +2 010 17 333 540
✉ info@nsbiotec.com



CMC Medical Devices &
Drugs S.L.
C/ Horacio Lengo, 18. 29006.
Malaga, Spain

General specifications:

- On-line monitoring to ensure good quality of water output
- Safe and reliable

Model	RO 100
Pretreatment cartridge	sediment filter 10 inch granular carbon filter 10 inch bacteria filter 10 inch A2 Solid carbon filter 10 inch
Membrane capacity	100 gallon/day
No. of pumps	1 pump
ultraviolet radiation unit	Less than 1 gallon/min
ion removal filter the risen tank	4*10
Reservoir tank	3 gallons
System display and indication parameters	guage for inlet pressure guage for outlet pressure T.D.S METER
Material	the cabinet made of treated iron sheet and coated with electrostatic paints (Stainless steel optional)
Outer Dimensions (W x D x H) mm	550 x 450 x 600
Power supply	220-240 VAC, 50Hz

